

PARASITOLOGY CENTER, INC. (PCI)

11445 E. Via Linda, # 2-419

Scottsdale, AZ 85259

Ph. (480) 767-2522

Fax (480) 767-5855

www.parasitetesting.com

INSTRUCTIONS FOR PINWORM TEST

Please read instructions completely before you begin.

1. Holding the cap, remove paddle from the container. Do not touch the paddle.
2. Separate the buttocks and press the sticky surface of the paddle against the anal area. The specimen should be obtained early in the morning before arising or bowel movement.
3. Return the paddle to the container and secure cap tightly.
4. Wash hands thoroughly with soap after you have finished.
5. Print patient information on the label.
6. Store specimen in a cool place until mailing back to PCI for analysis.
7. Results will be made to your health care practitioner within 2-5 working days.

If you have any questions please use the “contact us” link on our website www.parasitetesting.com

ONLY YOUR PHYSICIAN OR CLINIC IS AUTHORIZED TO DISCUSS YOUR RESULTS WITH YOU.